

EASTER BUSH SCIENCE OUTREACH CENTRE



**Get hands-on
with real-life
science**



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre

 www.ebsoc.ed.ac.uk
 @EBSOClab



A QUESTION OF TASTE

PCR Workshop



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre

Get hands-on
with real-life
science

Let's do a taste test!



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre





What are your taste test results?



Strong taster

- Quick negative reaction
- Very bitter
- *Ratings 4-5*



Weak taster

- Unsure at first, but then find the taste unpleasant
- *Ratings 1-3*



Non-taster

- Cannot taste anything bitter
- *Rating 0*


Record your own result on your sheet




THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre

Get hands-on
with real-life
science

What are you tasting?



PTC



Phenylthiocarbamide



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre



Get hands-on
with real-life
science

PTC is similar to the natural compounds found brussels sprouts



PTC only tastes bitter to around 70% of people. To the other 30% it is completely tasteless.

What was your phenotype?



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre

Get hands-on
with real-life
science

Why can only some of you taste it?



The ability to taste PTC has a **genetic basis**.



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre

Get hands-on
with real-life
science

What is your genotype?

Your **phenotype** is determined by your **genotype**.

Today, you will see if your genotype matches your phenotype.

Genotype = An organism's unique DNA

Phenotype = Observable characteristic of an organism



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre

Get hands-on
with real-life
science

What do we need to see your genotype?




THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre

Get hands-on
with real-life
science

Collecting your cheek cells



1. Write your lab number on the cup.
2. Swallow any excess saliva, clean your mouth with your tongue.
3. Gently chew the insides of your cheeks for 1 minute.



4. Swill your mouth with the salt water in the cup for 30 seconds. **Don't swallow it.**
5. Gently dribble liquid back into cup.

Make sure you only handle your own sample




Concentrate your cheek cells



1) Use marker pen to write your lab number on the **top and side** of a tube



2) Pipette the cell sample into the tube. Make sure you **take sample from bottom of cup and fill it to the top line (1.5ml)**



3) Close lid and place tube in small white foam rack, **one person** take it to the centrifuge.



Samples will be centrifuged to concentrate cells at bottom of tube




Heavier material (the cells) is thrown outwards and collects to form a **pellet**

Centrifuge spins tubes at 14,000 revolutions per minute (rpm) for 2 minutes


Liquid
(saliva, salt water)

Pellet
(cheek cells)



Get hands-on
with real-life
science

Discard the liquid



- 1) Pour most of the liquid **back into cup**.
- 2) Place the end of the tube on a **tissue** (to get out more liquid)
- 3) **Close lid**

Repeat if there are not enough cells.




Create a Cell Suspension (soup)



- 1) Break up the pellet of cells.
- 2) The cell suspension should be a **cloudy liquid with no lumps.**
- 3) When finished, place tube in your plastic rack.



Using Micropipettes- gather your equipment



P20 pipette

P200 pipette

Tips

Practice dye

Practice card



How to hold a micropipette



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre


Get hands-on
with real-life
science

Micropipette identification



1 - 20 μl

20-200 μl



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre

Get hands-on
with real-life
science

Setting the volume



1 - 20 µl

20-200 µl




Fig.1



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre




Let's practice setting the micropipettes!

$$1 \text{ ml} = 1000 \mu\text{l}$$



Can you set the pipette for 10 μl ?



Pipette set for 10 μl





Can you set the pipette for 30 μ l?

1 - 20 μ l 20-200 μ l

Pipette set for 30 μ l



Taking up liquid

- Put tip on pipette
- Press plunger down to 1st stop
- Place tip in liquid
- Release plunger SLOWLY



Dispensing liquid onto card

- Place tip where you want liquid to go
- Press plunger right down to 1st stop (when **putting liquid into a tube** we use the 2nd stop!)
- Move tip away from liquid
- Release plunger

Volume control
Plunger



Disposable
tip



Extract DNA from cells

1. Write your number on the lid **and** side of the screw cap tube containing Chelex beads.
2. Check that **micropipette** is set to 30 μ l.



3. Place a tip on pipette and cut the end of it at an angle.




30 μ l




4. Transfer 30 μ l of cell suspension **from** your sample tube **to** your new screw cap tube.

5. Screw lid down **tightly** on screw-cap tube, **flick** it and place in foam rack.

6. **One person** take it to the **vortex** and heat sample to break cell membranes and release DNA (100°C for 10 min)



What are the beads doing?




Chelex beads trap things which might
help cut the DNA




THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre

Get hands-on
with real-life
science

Chelex beads bind positive metal ions



So at the moment...



Chelex beads



heating




Centrifuge to remove used beads and cell debris



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre

Get hands-on
with real-life
science

Today we are looking at the gene TAS2R38



```
CCTTTCTGCACTGGGTGGCAACCAGGTCTTTAGATTAGCCAAGTAGAGAAGAGAAAGTAGA
ATAGCCAATTAGAGAAGTGACATCATGTTGACTCTAACTCGCATCCGCACTGTGTCCTAT
GAAGTCAGGAGTACATTTCTGTTTCATTTTCAGTCTGGAGTTTGCAGTGGGGTTTCTGACC
AATGCCTTCGTTTTTCTTGGTGAATTTTTGGGATGTAGTGAAGAGGCAGCCACTGAGCAAC
AGTGATTGTGTGCTGCTGTGTCTCAGCATCAGCCGGCTTTTCCTGCATGGACTGCTGTTT
CTCAGTTCCTAATCCAGCTTTAAGCAGTTTCAGTCAAGCAGTCAAGCAGCAGCTAC
```

This is the gene which encodes for a taste receptor protein that detects PTC


```
AACTCTCGTGACCCAGCCTGGAGGCCACATTAAGCCCTCAAGTCTCTTGTCTCCTTT
TTCTGCTTCTTTGTGATATCATCTGTGCTGCCTTCATCTCTGTGCCCTACTGATTCTG
TGGCGGACAAAATAGGGGTGATGGTTTGTGTTGGGATAATGGCAGCTTGTCCCTCTGGG
CATGCAGCCATCCTGATCTCAGGCAATGCCAAGTTGAGGAGAGCTGTGATGACCATTCTG
CTCTGGGCTCAGAGCAGCCTGAAGGTAAGAGCCGACCACAAGGCAGATTCCCGGACACTG
TGCTGAGAATGGACATGAAATGAGCTCTTCATTAATACGCCTGTGAGTCTTCATAAATAT
GCC
```

THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre




Does your tongue have taste receptors for PTC?


TAS2R38 taste receptor protein detects PTC




This is a single taste bud.



Today we're going to make billions of copies of this part of the gene



CCTTTCTGCACTGGGTGGCAACCAGGTCTTTAGATTAGCCAAGTATAGAGAAGAGAAGTAA
ATAGCCAATTAGAGAAGTGACATCATGTTGACTCTAACTCGCATCCGCACTGTGTCTTA
GAAGTCAGGAGTACATTTCTGTTTCAATTCAGTCCTGGAGTTTGCAGTGGGGTTTCTGACC
AATGCCCTTCGTTTTCTTGGTGAATTTTTGGGATGTAGTGAAGAGGCAGCCACTGAGCAAC



You will use **PCR** to make many copies of the part of the TAS2R38 gene we want to investigate.


CATGCAGCCATCCTGATCTCAGGCAATGCCAAGTTGAGGAGAGCTGTGATGACCATTCTG
CTCTGGGCTCAGAGCAGCCTGAAGGTAAGAGCCGACCACAAGGCAGATTCCCGGACACTG
TGCTGAGAATGGACATGAAATGAGCTCTTCATTAATACGCCTGTGAGTCTTCATAAATAT
GCC



The University of EDINBURGH
Easter Bush
Science Outreach Centre





What are the *ingredients* for PCR?



Setting up PCR reaction

Your will each receive a PCR tube containing



You will add the final ingredient:
your DNA!



Short break

Remove your lab coat if leaving the lab.




THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre



Setting up PCR reaction



1. Your tiny PCR tube is on ice




2. Write your number on lid **and** side, then place back on ice

3. Set pipette for 2.5 μ l






Add 2.5µl of your DNA to the PCR mix – check you have DNA!



Take DNA from **top of sample**, avoiding the beads







Centrifuge your tubes




Open the centrifuge



Put in the tubes, make sure they are balanced!



Close the centrifuge and turn on for 30 seconds then turn off



6 Tubes



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre

Get hands-on
with real-life
science

Your DNA now goes in the PCR machine



The PCR machine is programmed to run through all the required temperatures, and repeats the cycle 35 times



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre

Get hands-on
with real-life
science

A reminder about PCR



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre



Recreate the PCR reaction




G	A	C	C	G	A	T	G	A	C	A	G	A	T	G	A	T	T	A	G	G	A	A	A
C	T	G	G	C	T	A	C	T	G	T	C	T	A	C	T	A	A	T	C	C	T	T	T

DNA

C C G A T

A T C C T

primers



nucleotides



LUNCH



Please be back here for 1pm!



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre

Get hands-on
with real-life
science


What have we done so far?



Discovered our phenotype via PTC taste test



Sampled our own DNA



Burst open the cells and removed impurities from sample




Run a PCR to amplify a specific region of the PTC tasting gene



What will you do next?



Use molecular biology techniques to analyse your DNA




Compare the taste test results with your DNA analysis



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre

Get hands-on
with real-life
science

How does the TAS2R38 gene make us a....



Strong taster



Weak taster




Non-taster



A difference in just one nucleotide of *TAS2R38* gene

Point mutation at nucleotide 145



```
CCTTTCTGCACTGGGTGGCAACCAGGTCTTTAGATTAGCCAACCTAGAGAGAGAGAAGTAGA
ATAGCCAATTAGAGAAGTGACATCATGTTGACTCTAACTCGCATCCGCACTGTGTCTTAT
GAAGTCAGGAGTACATTTCTGTTTCAATTCAGTCCTGGAGTTTGCAGTGGGGTTTCTGACC
AATGCCTTCGTTTTCTTGGTGAATTTTTGGGATGTAGTGAAGAGGGGGCACTGAGCAAC
AGTGATTGTGTGCTGCTGTGTCTCAGCATCAGCCGGCTTTTCTGCAATGGACTGCTGTTC
CTGAGTGCTATCCAGCTTACCCACTTCCAGAAGTTGAGTGAACCACTGAACCACAGCTAC
CAAGCCATCATCATGCTATGGATGATTGCAAACCAAGCCAACCTCTGGCTTGCTGCCTGC
CTCAGCCTGCTTTACTGCTCCAAGCTCATCCGTTTTCTCTCACACCTTCTGATCTGCTTG
GCAAGCTGGGTCTCCAGGAAGATCTCCAGATGCTCCTGGGTATTATTCTTTGCTCCTGC
ATCTGCACTGTCCTCTGTGTTTGGTGCTTTTTTTAGCAGACCTCACTTCACAGTCACAAC
GTGCTATTCATGAATAACAATAACAAGGCTCAACTGGCAGATTAAAGATCTCAATTTATTT
TATTCCTTTCTCTTCTGCTATCTGTGGTCTGTGCCTCCTTTCTATTGTTTCTGGTTTCT
TCTGGGATGCTGACTGTCTCCCTGGGAAGGCACATGAGGACAATGAAGGTCTATAACCAGA
AACTCTCGTGACCCAGCCTGGAGGCCACATTAAAGCCCTCAAGTCTCTTGTCTCCTTT
TTCTGCTTCTTTGTGATATCATCCTGTGCTGCCTTCATCTCTGTGCCCTACTGATTCTG
TGGCGGACAAAATAGGGGTGATGGTTTGTGTTGGGATAATGGCAGCTTGTCCCTCTGGG
CATGCAGCCATCCTGATCTCAGGCAATGCCAAGTTGAGGAGAGCTGTGATGACCATTCTG
CTCTGGGCTCAGAGCAGCCTGAAGGTAAGAGCCGACCACAAGGCAGATTCCCGGACACTG
TGCTGAGAATGGACATGAAATGAGCTCTTCATTAATACGCCTGTGAGTCTTCATAAATAT
```



UNIVERSITY of EDINBURGH

Edinb

Science Outreach Centre

Get hands-on
with real-life
science

How many chromosomes do you have?

What do we call
different forms
of the same
gene?



Alleles


Humans have 23 pairs of chromosomes
→ We have **2 copies of each gene**





THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre

Get hands-on
with real-life
science

Alleles are different forms of the same gene



Same alleles:




Different alleles:



Alleles can be dominant or recessive

 Allele for not tasting PTC (t)
 Allele for tasting PTC (T)

Make a prediction...
what do you think
your genotype is?



PTC taster
allele is
dominant

GENOTYPE

Tt
Heterozygous

TT
Homozygous

tt
Homozygous

PHENOTYPE



Weak taster

Strong taster

Non-taster



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre

Get hands-on
with real-life
science

How can we tell if our DNA contains a C or G at position 145?

A difference in just one nucleotide of TAS2R38 gene

Point mutation at nucleotide 145

CCITTTCTGCACTGGGTGGCAACCAGGTCTTTAGATTAGCCAAGTACAGAGAGAAGTAGA
ATAGCCAATTAGAGAAGTGACATCATGTTGACTCTAACTCGCATCCGGCTGTGTCCTAT
GAAGTCAGGAGTACATTCTGTTCATTTGAGTCTGAGTTTGCAGTGGGCTTCTGACC
AATGCCTTCGTTTTCTTGGTGAATTTTGGGATGTAGTGAAGAGGGCCCACTGAGCAAC
AGTGATTGTGTGCTGCTGTCTCAGCATCAGCCGGCTTTTCTGCTGGACTGCTGTTC
CTGAGTGCTATCCAGCTTACCCACTTCCAGAAGTTGAGTGAACCACTGAACCACAGCTAC
CAAGCCATCATCATGCTATGGATGATTGCAAACCAAGCCAACCTTGGCTTGCTGCCTGC
CTCAGCCTGCTTTACTGCTCCAAGTCTCCAGATGCTCCTGGGTATTATTCTTGGCTCCTGC
GCAAGCTGGGTCTCCAGGAAGATCTCCAGATGCTCCTGGGTATTATTCTTGGCTCCTGC
ATCTGCACTGTCTCTGTGTTTGGTCTTTTAGCAGACTCACTTCCAGTCAACA
GTGCTATTTCATGAATAACAATAAAGGCTCAACTGGCAGATTAAAGATCTCAATTTATT
TATTCCTTCTCTCTGCTATCTGTGGTCTGTGCTCCTTCTTCTGTTTCTGGTTTCT
TCTGGGATGCTGACTGTCTCCCTGGGAAGGCACATGAGGACAATGAAGTCTATAACCAGA
AACTCTCGTGACCCAGCCTGGAGGCCACATTAAAGCCCTCAAGTCTCTGTCTCCTTT
TTCTGCTTCTTGTGATATCATCTGTGCTGCCTCATCTCTGTGCCCTACTGATTCTG
TGGCCGACAAAATAGGGGTGATGTTGTTGGGATAAAGGAGCTGTGATGACCATTCTG
CATGCAGCCATCCTGATCTCAGGCAATGCCAAGTTGAGGAGAGCTGTGATGACCATTCTG
CTCTGGGCTCAGAGCAGCTGAAGGTAAGAGCCGACCACAAGGCAGATCCCGGACACTG
TGCTGAGAAATGGACATGAAATGAGCTCTTCAATAACGCCTGTGAGTCTTCATAAATAT

Get hands-on with real-life science

A difference in just one nucleotide of TAS2R38 gene

Point mutation at nucleotide 145

CCTTTCTGCACTGGGTGGCAACCAGGTCTTTAGATTAGCCAAGTACAGAGAGAAGTAGA
ATAGCCAATTAGAGAAGTGACATCATGTTGACTCTAACTCGCATCCGGCTGTGTCCTAT
GAAGTCAGGAGTACATTCTGTTCATTTGAGTCTGAGTTTGCAGTGGGCTTCTGACC
AATGCCTTCGTTTTCTTGGTGAATTTTGGGATGTAGTGAAGAGGGCCCACTGAGCAAC
AGTGATTGTGTGCTGCTGTCTCAGCATCAGCCGGCTTTTCTGCTGGACTGCTGTTC
CTGAGTGCTATCCAGCTTACCCACTTCCAGAAGTTGAGTGAACCACTGAACCACAGCTAC
CAAGCCATCATCATGCTATGGATGATTGCAAACCAAGCCAACCTTGGCTTGCTGCCTGC
CTCAGCCTGCTTTACTGCTCCAAGTCTCCAGATGCTCCTGGGTATTATTCTTGGCTCCTGC
GCAAGCTGGGTCTCCAGGAAGATCTCCAGATGCTCCTGGGTATTATTCTTGGCTCCTGC
ATCTGCACTGTCTCTGTGTTTGGTCTTTTAGCAGACTCACTTCCAGTCAACA
GTGCTATTTCATGAATAACAATAAAGGCTCAACTGGCAGATTAAAGATCTCAATTTATT
TATTCCTTCTCTCTGCTATCTGTGGTCTGTGCTCCTTCTTCTGTTTCTGGTTTCT
TCTGGGATGCTGACTGTCTCCCTGGGAAGGCACATGAGGACAATGAAGTCTATAACCAGA
AACTCTCGTGACCCAGCCTGGAGGCCACATTAAAGCCCTCAAGTCTCTGTCTCCTTT
TTCTGCTTCTTGTGATATCATCTGTGCTGCCTCATCTCTGTGCCCTACTGATTCTG
TGGCCGACAAAATAGGGGTGATGTTGTTGGGATAAAGGAGCTGTGATGACCATTCTG
CATGCAGCCATCCTGATCTCAGGCAATGCCAAGTTGAGGAGAGCTGTGATGACCATTCTG
CTCTGGGCTCAGAGCAGCTGAAGGTAAGAGCCGACCACAAGGCAGATCCCGGACACTG
TGCTGAGAAATGGACATGAAATGAGCTCTTCAATAACGCCTGTGAGTCTTCATAAATAT


Get hands-on with real-life science



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre

Get hands-on
with real-life
science

We can use restriction enzymes



Restriction enzymes cut the DNA at specific sequences



Different restriction enzymes cut different sequences



EcoRI

G A A T T C
C T T A A G

PstI

C T G C A G
G A C G T C




G A A T T C A C G T C T G C A G C C A A A T G G C G A A T T C C A
C T T A A G T G C A G A C G T C G G T T T A C C G C T T A A G G T



Which allele with HaeIII cut?


GGCC



HaeIII


Taster

ggggggggggccaatgggga



Non taster

gaggcgggcaactgagca



Safety first!




THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre



Setting up the enzyme reaction



1. Add 10 μ l of your PCR product to bottom of each tube



2. To **W**, Add 5 μ l **water**

*Change pipette tips
between tubes!*


3. To **E**, Add 5 μ l **enzyme**

4. Flick the tubes then centrifuge for about 30 sec and put it in a foam rack.

5. **One person** take the tubes to the water bath and then incubate at 37°C for 30 minutes




How can we tell if the enzyme has cut the DNA?




Gel Electrophoresis

Gel electrophoresis



Group activity



GGCC
CCGG





HaeIII

Loading sample on to gel



Hold pipette tip
just above well,
below buffer
level




Be careful not to pierce
bottom of well with pipette tip!

This time you push down to the 1st stop to fill the well with the DNA

Have a go!



Load 10 μ l practice dye (tube marked **P**) to each well.



Push down to the 1st stop to fill the well with the DNA



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre

Get hands-on
with real-life
science

Prepare the DNA for analysis




1. Add 2.5 μ l orange loading dye (tube marked LD) to sample W and E
2. **Change the pipette tip** for each sample
3. Flick then centrifuge samples (**make sure they are balanced**)



Meet the Scientists




Loading DNA on to gel



Which allele with HaeIII cut?


GGCC



HaeIII


Taster

ggggggggggccaatgggga




Non taster

gagggcgggcactgagca



Possible digest results



Practice your analysis...



1) What is the phenotype of the tongue?

Taster

2) What is the genotype of the tongue?

TT or Tt

3) Will the enzyme cut?

Yes


4) How many fragments of DNA will there be?

If TT then two fragments

If Tt then three fragments



This is what it will really look like...



Safety check!



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre






Viewing the DNA



A DNA stain was added to the agarose when the gels were made



Unexpected results?

Phenotype	Expected Genotype	Actual Genotype	Possible reasons
	TT	Tt	<ul style="list-style-type: none">• High density of taste buds• Other genes involved
	Tt	TT	<ul style="list-style-type: none">• Low density of taste buds• Dry mouth• Have a cold• Other genes involved
	tt	Tt	<ul style="list-style-type: none">• Low density of taste buds• Dry mouth• Have a cold• Other genes involved



EASTER BUSH SCIENCE OUTREACH CENTRE



**Get hands-on
with real-life
science**




THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre

 www.ebsoc.ed.ac.uk
 @EBSOClab

Get hands-on
with real-life
science

Please write down three
words on our wall that
describe your experience
today!



THE UNIVERSITY of EDINBURGH
Easter Bush
Science Outreach Centre



Fun

Boring

Informative

Inspiring

Rewarding

Uninteresting

Interesting

Confusing

Enjoyable

Difficult

Thought-provoking

Frustrating

Dull

